



VITAL

Connecting providers and patients around the globe™



Blood Pressure Monitor

Convenient blood pressure and pulse readings

DESCRIPTION

The VITAL Blood Pressure Monitor is a **clinically validated** device that utilizes the IEEE-standard Bluetooth communications protocol to capture and transmit systolic and diastolic blood pressure readings as well as pulse readings to the VITAL Tablet using a wireless connection.

The blood pressure monitor also includes voice prompts that guide the patient through each step involved in taking a blood pressure reading, including placing the cuff on the arm, pressing the "Start" button to take the blood pressure reading, and transmitting the reading automatically to the VITAL tablet.

Once the reading is captured on the VITAL Tablet, then the reading may be shared with the clinician through the VITAL Medical Video (MeVi) application.

KEY FEATURES



Medium Cuff: 9.0 - 14.6 inch for gentle inflation



Automatic irregular heartbeat detector



Bluetooth connectivity



Easy-to-read LCD screen



Power source: 4 AAA batteries



Measurement range
Pressure: 20-280 mmHg
Systolic pressure: 60-279 mmHg
Diastolic pressure: 40-200 mmHg
Pulse: 40-180 beats/min



Measurement accuracy
Pressure: 3 mmHg
Pulse: 5%

Contact us today to find out how Blood Pressure Monitor can benefit your organization!